

Gryllus Techniques and Protocols**Total RNA Extraction & cDNA Synthesis**

To extract RNA from embryos and make a library of cDNA from cricket embryos, as a more stable alternative to mRNA.

General Notes:

- Be certain to clean bench and pipettes first with water and 70% ethanol.
- Use ~20 embryos per tube.
- Can repeat cDNA reaction with 2x, 3x, 4x etc volume if the reaction fails.

- 1) Add 200µl Trizol.
*Note: Steps involving Trizol must be performed under a fume hood and with caution, as Trizol is toxic.
- 2) Homogenize with pestle, by hand.
- 3) Add 800µl Trizol; rinse pestle over tube as you pipette, to collect sample stuck to pestle.
- 4) Spin 10min at 12,000rpm at 4°C in a microfuge.
- 5) Incubate 5min at room temperature.
- 6) Add 200µl chloroform.
* Note: Steps involving chloroform must be performed under a fume hood and with caution, as chloroform is toxic.
- 7) Vortex for 15sec.
- 8) Incubate 3min at room temperature.
- 9) Spin 15min at 12,000rpm at 4°C in a microfuge.
- 10) Collect supernatant and add 500µl isopropanol.
- 11) Spin 15min at 12,000rpm at 4°C in a microfuge.
- 12) Aspirate and discard s/n.
- 13) Wash with 500µl 70% EtOH in MilliQ H₂O.
- 14) Spin 15min at 12,000rpm at 4°C in a microfuge.
- 15) Aspirate and discard s/n.
- 16) Air dry 2-3min. Do not fully dry, as it can be difficult to dissolve completely dried RNA.
- 17) Resuspend in 10µl MilliQ H₂O.
Use water from a new tube, that is certainly clean.
- 18) Measure concentration with spectrophotometer.
- 19) Use Invitrogen SuperScript III First-Strand Synthesis System for RT-PCR to generate cDNA.
- 20) Measure concentration with spectrophotometer. Can also perform a PCR and gel electrophoresis with reliable primers/conditions to test cDNA quality.